

# OxyHOME<sup>®</sup>

## Stationary Concentrator

Keep Going.<sup>SM</sup>

NEW

### BENEFITS

- Quiet Operation
- Small Footprint
- Purity Check by Percentage
- Alarm Indicators: No Flow, Low/High Pressure, Purity, Power
- German Compressor
- Twin Blower Cooling Design
- Contemporary Design
- 3-Year Replacement Warranty

### PRODUCT SPECS

Part Number 1400-5100

Flow Specifications	0 - 5 LPM
Oxygen Concentration	93% (+3%, -3%)
Oxygen Output Pressure	8 (+2%, -2%) psi
Rated Voltage/Frequency	120 VAC ±10% 60Hz
Power Consumption	350 Watts
Sound	≤ 45 dB(A)
Dimensions	13.4" W x 11.8" D x 25.6" H 34 cm x 30 cm x 65 cm
Weight	41.8 lbs (19 kg)



# 1400-5100

CONTACT US TODAY TO LEARN MORE

888-327-7301 · [orders@oxygo.life](mailto:orders@oxygo.life) · [oxygo.life](http://oxygo.life)